Hemospon[®]

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Hemospon[®]

Hemostatic Sponge made with
Lyophilized Hydrolyzed Collagen (Gelatin)

Indication: Hemospon is indicated during and after surgical procedures for achieving local homeostasis, especially when the control of bleeding by conventional methods is not efficient. During dental surgical procedures, Hemospon, besides allowing hemostasis, also has the usages listed below:

- Hemospon 1x1x1 cm: acts in the maintenance of blood clot and in the filling of the spaces generated en tooth extractions, removal of impacted tooth, removal of cysts or tumors, biopsies, etc.
- Hemospon Tape e Cote: protection of wounds produced during dental procedures, such as periodontal surgeries, closure of graft areas and maxillary sinus surgery.
- Hemospon Size and Standard: Can be applied in general surgery, plastic surgery, dermatology, gynecology, dentistry, otolaryngology, proctology, etc

Composition: Gelatin of pig origin

Description of the Product: Hemospon is a lyophilized gelatin. Hemospon has hemostatic and healing action and it is completely absorbed by the organism within approximately 15 days. Hemospon is non-toxic and non-pirogenic. It is presented ready for use, packaged in individual blisters. The product is presented sterile by irradiation.

Mechanism of action: Hemospon should be applied as an absorbable hemostatic, used to reduce bleeding and stabilize the blood clot. Its porous structure allows the product to have a high capacity of fluid absorption. It absorbs 40 to 50 times its weight of liquid. It is well known that blood clot forms due to the platelet's agglutination in the place of the bleeding and also due to the activation of the coagulation factors cascade. Hemospon works as a matrix that allows additional reinforcement to the formed blood clot. Besides this local effect, the product has no other pharmacodynamic effects.

Instructions for use:

Remove Hemospon from packaging, using aseptic technique. If required, the sponge may be cut to a desired size.

Apply Hemospon directly to the area where hemostasis is required, making a slight pressure.

The product may be used dry, or in option, moistened in sterile saline solution.

The necessary time for achieving hemostasis depends on the intensity of bleeding.

Once hemostasis is achieved, Hemospon may be removed, substituted or left in place until its complete absorption, which happens within 15 days.

Contra-indications: The usage of sponges is not recommended in infected or contaminated areas.

Warning and special care: Hemospon is presents sterile and should not be re-sterilized. Before usage, check the package. Do not use if package is damaged. Single use product. The units must be handled with sterilized tweezers or scissors. There is no report f adverse effect for the product used correctly according to its indications.

Storage: Keep at room temperature

Sterile product – Sterilized by irradiation.

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PRODUTOS ODONTOLÓGICOS S.A.